

**Medical Cool Box Product Number: MCB**

<p><b>High capacity fast cooling intelligent inverter-controlled system. This professional portable medical cooler can store between 3 and 5 insulin injector pens or 9 insulin or about 20 refills.</b></p> <p><b>It keeps your medicine at a cool 2° - 8° Celsius in 30°C ambient environment with a temperature range of between 2° - 25° Celsius.</b></p> <p><b>A long 12 hour stand-by is provided from its powerful built-in 10400mAh Samsung battery. Ready for travel with the provided car charger and portable bag inclusive when you purchase this outstanding medical cooler.</b></p>	
---	--

<b>Internal Dimensions:</b>	175mm x 60mm x 26mm
<b>External Dimensions:</b>	207mm x 94mm x 91mm
<b>Battery</b>	<b>Built-in 10400mAh Samsung battery.</b>
<b>Battery Life</b>	Internal battery = 12 Hour Stand by External battery = 10 Hours stand by TOTAL = 22 Hours
<b>Included</b>	Car charger, portable bag, 1 X internal battery, 1x external battery, USB charge cable, instruction manual 1 x UK adaptor, 1 x European adaptor
<b>Weight (KG)</b>	5 KG