

PRODUCT DATA SHEET JYK-CIB

MEDICAL VACCINE COOLER: Product Number: JYK-CIB

Description:

This cool box is perfect for carrying Personal Medication, Insulin, Vaccines, Botox Vials, Blood samples, Biological samples and heat-sensitive medication distribution. It is also suitable for Health Centres and Animal Husbandry. CE Certified.

This product can store 3-4 injection pens (Novo pen) or 3 one-time injection pens; Insulin and interferon can be hold about 9 teams or 15-20 leads.



(Pens not included)

TEMPERATURE	Keeps 2-8 degrees C under ambient 30 degs C
Weight Empty:	0.65kg
Capacity/Internal Dimensions:	175 x 60 x 25mm
External Dimensions:	210 x 90 x 95 mm
Refridgeration Temperature	2-8 °C (under ambient 320C)
Input Voltage	DC7.4-12V
Product Power	8W-24W +- 5%
Alarm	Yes
Tunable Temperature	2-25 C
Battery	7.4V 4000mAh Lithium
Display	LED temperature display, easy button operation
Accessories	Battery, domestic adaptor, power cord, instructions, warranty card.
Warranty	One year

The Cool Ice Box Company Ltd.
 27 Broad St, South Molton , Devon EX36 3AQ
 UK Office no: 01598 740 685
 European Office: 00 33 (0) 553836230

www.coolicebox.com
info@coolicebox.com
**NOMAD is the TM of The Cool Ice
 Box Company Ltd.**

Registered Address:
 13 Mill Street
 Bideford, Devon
 EX39 2JJ