

PRODUCT DATA SHEET SB1215

Medical Cold Chain Cooler Bag Product Series: Product Number: SB1215

Description:

SB1215 is a light weight Medical Cold Chain Soft Sided Cooler Bag designed specifically to be used principally as a stand alone temporary containment or storage function where safe temperature limits are essential in transit, or during use, for medical products including vaccines for human beings in the prevention, treatment or alleviation of disease. SB1215 is in conformity with the Medical Device Directive 93/42/EEC for CE Marking.



Usable Cold Space Capacity/Dimensions:	8.1L (2.1 US Gal.)/25x17x19cm (9.8x6.7x7.5")
Weight Empty/Weight Incl. Gel packs:	0.925kg (2.0 lb.)/3.32kg (7.3 lb.)
Empty Capacity/Internal Dimensions:	12.3L (3.3 US Gal.)/28x20x22cm (11x7.9x8.7")
External Dimensions:	30x22x24cm (11.8x8.7x9.5")
External Surface/Internal Surface:	Polyester 600D + EVA thermoform sheet front & bottom/PVC Waterproof
Liner:	PVC transparent waterproof with Gel pockets
Insulation:	Sides & Bottom 10mm (0.39") NBR Closed cell foam (K=0.034)
Lid Type:	Polyester 600D. Insulated, lined, fabric hinge, 2 x heavy duty Zips
Refrigerant Pack/No. of/Dimensions:	HTI21 Gel (reusable)/12/21x8cm (8.3x3.2"), 200g (7.1oz.) +/-10%
Digital Thermometer:	Optional
Temperature Data Recorder:	Optional
Packing:	PE bag + cardboard outer carton
Shipping Carton Dimensions: Pieces/Weight per Carton:	48x36x31.5cm (18.9x18.1x12.4") 4/14.5kg (32.0lbs)
Cartons/Pcs. per 20' TEU (Total pcs.):Net Weight (FCL):	400/1600 5,800kg (6.4 US Short Tons)
Cold Life from Close of Lid to +10oC:	18 hrs 20 mins. @ 30oC (86oF) Ambient (User may specify Cold life requirement. Usable Cold Space may change accordingly.)
Minimum Temperature in Cold Space:	+2oC (35.6oF)

NOTES:
 UV resistant external material;Product does not contain Lead, Mercury, Cadmium, Hexavalent Chromium,Polybrominated Biphenyls (PBB), Polybrominated Biphenyl Ethers (PBBE) or compounds containing Benzene
 The closed cell foam insulation used does not use CFC as the blowing agent
 Ingress protection meets IEC 60529 - IP55

The Cool Ice Box Company Ltd.
 27 Broad St, South Molton , Devon EX36 3AQ
 UK Office no: 01598 740 685
 European Office: 00 33 (0) 553836230
 Mobiles: 07968 797785 07710 471598
VAT No. 118 9513 02

www.coolicebox.com
info@coolicebox.com
**NOMAD is the TM of The Cool Ice
 Box Company Ltd.**

Registered Address:
 13 Mill Street
 Bideford, Devon
 EX39 2JJ