

**PRODUCT DATA SHEET  
SB1215T**

**Medical Cold Chain Cooler Bag Product Series: Product Number: SB1215T**

**Description:**

SB1215T is a light weight Medical Cold Chain Soft Sided Cooler Bag with fitted LED display showing the internal temperature - designed specifically to be used principally as a stand alone temporary containment or storage function where safe temperature limits are essential in transit, or during use, for medical products including vaccines for human beings in the prevention, treatment or alleviation of disease.

SB1215T is in conformity with the Medical Device Directive 93/42/EEC for CE Marking.



<b>Usable Cold Space Capacity/Dimensions:</b>	8.1L (2.1 US Gal.)/25x17x19cm (9.8x6.7x7.5")
<b>Weight Empty/Weight Incl. Gel packs:</b>	0.925kg (2.0 lb.)/3.32kg (7.3 lb.)
<b>Empty Capacity/Internal Dimensions:</b>	12.3L (3.3 US Gal.)/28x20x22cm (11x7.9x8.7")
<b>External Dimensions:</b>	30x22x24cm (11.8x8.7x9.5")
<b>External Surface/Internal Surface:</b>	Polyester 600D + EVA thermoform sheet front & bottom/PVC Waterproof
<b>Liner:</b>	PVC transparent waterproof with Gel pockets
<b>Insulation:</b>	Sides & Bottom 10mm (0.39") NBR Closed cell foam (K=0.034)
<b>Lid Type:</b>	Polyester 600D. Insulated, lined, fabric hinge, 2 x heavy duty Zips
<b>Refrigerant Pack/No. of/Dimensions:</b>	HTI21 Gel (reusable)/12/21x8cm (8.3x3.2"), 200g (7.1oz.) +/- 10%
<b>Digital Thermometer:</b>	Included
<b>Temperature Data Recorder:</b>	Optional
<b>Packing:</b>	PE bag + cardboard outer carton
<b>Shipping Carton Dimensions: Pieces/Weight per Carton:</b>	48x36x31.5cm (18.9x18.1x12.4") 4/14.5kg (32.0lbs)
<b>Cartons/Pcs. per 20' TEU (Total pcs.):Net Weight (FCL):</b>	400/1600 5,800kg (6.4 US Short Tons)
<b>Cold Life from Close of Lid to +10oC:</b>	18 hrs 20 mins. @ 30oC (86oF) Ambient (User may specify Cold life requirement. Usable Cold Space may change accordingly.)
<b>Minimum Temperature in Cold Space:</b>	+2oC (35.6oF)

**NOTES:**

UV resistant external material;Product does not contain Lead, Mercury, Cadmium, Hexavalent Chromium, Polybrominated Biphenyls (PBB), Polybrominated Biphenyl Ethers (PBBE) or compounds containing Benzene  
The closed cell foam insulation used does not use CFC as the blowing agent

**The Cool Ice Box Company Ltd.**  
27 Broad St, South Molton, Devon EX36 3AQ  
UK Office no: 01598 740 685  
European Office: 00 33 (0) 553836230  
Mobiles: 07968 797785 07710 471598  
VAT No. 118 9513 02

[www.coolicebox.com](http://www.coolicebox.com)  
[info@coolicebox.com](mailto:info@coolicebox.com)  
**NOMAD is the TM of The Cool  
Ice Box Company Ltd.**

Registered Address:  
13 Mill Street  
Bideford, Devon  
EX39 2JJ