

## PRODUCT DATA SHEET SB7003

### Nomad Medical Cold Chain Cooler Product Number SB7003

<p><b>Description:</b> SB7003 is a soft sided Cold Chain Cooler designed for low temperature storage of Vaccines and Medical Products, including IVD and Biological Specimens.</p> <p><b>Manufacturer's Warranty:</b> The Manufacturer warrants all products free from defects in material and workmanship under normal use and service for a period of one year from the date of purchase. This warranty shall not apply to any product that has been (i) subject to misuse, abuse or neglect; (ii) damaged accidentally; (iii) used for purposes not expressly stated; or (iv) repaired or altered in any way. Warranty Claims must be made prior to the expiry date of the warranty period.</p>	
--	--

<b>Usable Cold Space Capacity/Dimensions:</b>	12L (3.17 US Gal)/26x19x24 cm. (10.23x7.48x9.45 inches)
<b>Weight Empty/Weight Incl. Gel packs:</b>	0.5 kg. (1.1 lb.)/4.1 kg. (8.6lb.)
<b>Empty Capacity/Internal Dimensions:</b>	18L (4.76 US Gal)/L28xD23xH28 cm. (11.02x9.06x11.02 inches)
<b>External Dimensions:</b>	L30xD25xH30 cm. (11.81x9.84x11.81 inches)
<b>External Surface/Internal Surface:</b>	Polyester 600D/PVC (waterproof)
<b>Liner:</b>	Nylon 420D (removable)
<b>Lid Type:</b>	Polyester 600D, Insulated, Lined, Fabric Hinge, Heavy Duty Zip
<b>Insulation:</b>	8 mm. (0.32 inch) PE closed cell foam
<b>Refrigerant Pack/No. of/Dimensions:</b>	HTIceGel (Re-usable)/28x8 cm. (11.02x3.15 inches)
<b>Number/Mass of Refrigeration Packs:</b>	12/300 g.+or-20 g.(10.72 oz.+or-0.71 oz.)
<b>Digital Thermometer:</b>	Optional
<b>Temperature Data Recorder:</b>	Optional
<b>Cold Life from Close of Lid to 8oC (46.4oF):</b>	18 hrs. 21 mins.
<b>Cold Life from Close of Lid tp 10oC (50.0oF):</b>	19 hrs.
<b>Packing:</b>	PE bag + Cardboard outer carton
<b>Shipping Carton Dimensions:</b>	46.5x40.25x36.5 cm.(18.31x15.85x14.37 inches)
<b>Pieces/Weight per Carton:</b>	4/17.8 kg. (39.24 lbs.)
<b>Compliance:</b>	CE EC Medical Devices Directive 93/42/EEC, Class I

[www.coolicebox.com](http://www.coolicebox.com)

[info@coolicebox.com](mailto:info@coolicebox.com)

**NOMAD is the TM of The Cool Ice Box Company Ltd.**

**The Cool Ice Box Company Ltd.**

27 Broad St, South Molton , Devon EX36 3AQ

UK Office no: 01598 740 685

European Office: 00 33 (0) 553836230

Mobiles: 07968 797785 07710 471598

Registered Address:  
13 Mill Street  
Bideford, Devon  
EX39 2JJ