

PRODUCT DATA SHEET TCSS

Nomad Medical Cold Chain Cooler Product Number TCSS

<p>TCSS are low temperature holding packs for Insulin and designed specifically by the manufacturer to be used principally as a stand alone temporary containment or storage function within safe temperature limits for the contents while in transit or during use, for:</p> <ol style="list-style-type: none"> 1. Insulin for human beings for the purpose of the treatment or alleviation of Diabetes Type 1 and 2; and 2. Syringes, disposable Insulin Pen Systems and Blood Glucose Monitor. 	
--	--

Cold Life at 35o C ambient:	3 to 5 hrs. (2 -10oC)
Example of Contents (Not Supplied)	Insulin - 3*10 ml vials Disposable Insulin Pen Systems - 2 Syringes - 3 Blood Glucose Monitor - 1
Usable Cold Space Capacity/Dimensions:	190x100x35 mm
Weight Empty/Weight Incl. Gel packs:	130 gm/230g
Storage Method:	Net Pocket (zipped), Elastic webbing storage loops
External Dimensions:	200x110x45 mm
External Surface:	Ethylene Vinyl Acetate (EVA)
Internal Surface:	420 Denier Nylon
Lid Type:	Sealed hinge and zip fastener, with inner lip to reduce cold losses
Insulation:	PE closed cell foam - 5.0 mm (k = 0.076 mW/m-K @ 0oC)
Refrigerant Pack/No. of/:	1 off HTIceGel Pack, HTI-19
HTIceGel Pack Dimensions:	190x80x10 mm
HTIceGel Pack Weight:	100 gms/pc
Digital Thermometer:	Recommended Optional
Shipping Carton Dimensions: Pieces/Weight per Carton:	42x42x58cm 100pc

The Cool Ice Box Company Ltd.

27 Broad St, South Molton , Devon EX36 3AQ

UK Office no: 01598 740 685

European Office: 00 33 (0) 553836230

Mobiles: 07968 797785 07710 471598

www.coolicebox.com

info@coolicebox.com

**NOMAD is the TM of The Cool Ice
Box Company Ltd.**

Registered Address:
13 Mill Street
Bideford, Devon
EX39 2JJ